

## INDIAN INSTITUTE OF TECHNOLOGY PATNA DEPARTMENT OF CHEMISTRY

Bihta, Amhara Road, Patna – 801106, India; Phone +91-612 3028225

Sample submission form

To The HoD(Chemistry Department) IIT Patna, Bihta Amhara Road Patna-801106

**Subject: Usage of SDT Simultaneous TGA (Q600)** 

Name	:		
Designation	:		
Address for communica	tion :		
Phone number	:		
Email Address	:		
DD No	Date	Amount	Bank
~			

Certified that the sample submitted belong to the user mentioned above.

Signature of Guide with date & Seal Signature with date & Seal

(HOD/Principal/Managing director)

Date:

## **Sample information**

## Sample ID

Programme steps for heating/cooling Temperature Range

**Composition or nature of samples** 

Is there any expected reaction with the crucible or melting of any component in the temperature range of analysis

Atmosphere TGA Nitrogen

Sample pan Alumina cup 90µl

Charges for the measurement should be sent through an advance Demand Draft in favour of 'IIT Patna' payable at 'Patna' along with the sample or through RTGS/NEFT( Please mention A/c holder name......, A/c No......,Bank name.....,IFSC code......) to HOD, Chemistry Dept., IIT Patna (Please attach a copy of the transaction details along with signature of guide with date and seal)